

Product Datasheet

PPIL3 antibody (orb140367)

Description	Rabbit Polyclonal to PPIL3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PPIL3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: PPIL3. Gene-IDs: 53938. Accessions: NP_115861.1, NP_570981.1, Q9H2H8.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18155
Uniprot ID	Q9H2H8
NCBI	NP_570981.1, NP_115861.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

