



Product Datasheet

EYA2 antibody (orb140383)

Description	Rabbit Polyclonal to EYA2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	EYA2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: EYA2. Gene-IDs: 2139. Accessions: O00167, NP_742108.2, NP_005235.3.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59232
Uniprot ID	O00167
NCBI	NP_005235.3 , NP_742108.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...