

Product Datasheet

POLR1E antibody (orb140388)

Description	Rabbit Polyclonal to POLR1E.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	POLR1E
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: POLR1E. Gene-IDs: 64425. Accessions: NP_071935.1, Q9GZS1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47260
Uniprot ID	Q9GZS1
NCBI	NP_071935.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

