

Product Datasheet

Prolactin antibody (orb140393)

Description	Rabbit polyclonal antibody to Prolactin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PRLHR
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: PRLHR. Gene-IDs: 2834. Accessions: P49683, NP_004239.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41121
Uniprot ID	P49683
NCBI	NP_004239.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

