

Product Datasheet

GFOD2 antibody (orb140412)

Description	Rabbit Polyclonal to GFOD2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GFOD2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: GFOD2. Gene-IDs: 81577. Accessions: NP_110446.3, Q3B7J2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42255
Uniprot ID	Q3B7J2
NCBI	NP_110446.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

