

Product Datasheet

KLHL3 antibody (orb140415)

Description	Rabbit Polyclonal to KLHL3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KLHL3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: KLHL3. Gene-IDs: 26249. Accessions: Q9UH77, NP_059111.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64970
Uniprot ID	Q9UH77
NCBI	NP_001244124.1, NP_059111.2, NP_001244123.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

