

Product Datasheet

SYT17 antibody (orb140418)

Description	Rabbit Polyclonal to SYT17.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SYT17
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: SYT17. Gene-IDs: 51760. Accessions: Q9BSW7, NP_057608.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53849
Uniprot ID	Q9BSW7
NCBI	NP_057608.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

