

Product Datasheet

CDS2 antibody (orb140423)

Description	Rabbit Polyclonal to CDS2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CDS2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: CDS2. Gene-IDs: 8760. Accessions: O95674, NP_003809.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51418
Uniprot ID	O95674
NCBI	NP_003809.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

