

Product Datasheet

CLP1 antibody (orb140425)

Description	Rabbit Polyclonal to CLP1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CLP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: CLP1. Gene-IDs: 10978. Accessions: NP_006822.1, Q92989, NP_001136069.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47646
Uniprot ID	Q92989
NCBI	NP_006822.1 , NP_001136069.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

