

Product Datasheet

CRYZL1 antibody (orb140427)

Description	Rabbit Polyclonal to CRYZL1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CRYZL1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: CRYZL1. Gene-IDs: 9946. Accessions: O95825, NP_665857.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38697
Uniprot ID	O95825
NCBI	NP_665857.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

