

Product Datasheet

PAFAH1B2 (Phospho-Ser64) antibody (orb140435)

Description	Rabbit polyclonal to PAFAH1B2 (Phospho-Ser64)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	PAFAH1B2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: PAFAH1B2. Gene-IDs: 5049. Accessions: NP_001171677.1, NP_002563.1, P68402.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25569
Uniprot ID	P68402
NCBI	NP_001171676.1, NP_001171675.1, NP_002563.1, NP_001171677.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

