

Product Datasheet

MS4A1 antibody (orb140444)

Description	Mouse Monoclonal to MS4A1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	MS4A1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: MS4A1. Gene-IDs: 931. Accessions: NP_690605.1, P11836, NP_068769.2.
Isotype	IgG1
Clonality	Monoclonal
MW	33077
Uniprot ID	P11836
NCBI	NP_068769.2, NP_690605.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

