

Product Datasheet

LGR5 antibody (orb140452)

Description	Mouse Monoclonal to LGR5.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	LGR5
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Application notes	Gene Symbol: LGR5. Gene-IDs: 8549. Accessions: NP_003658.1, O75473.
Isotype	IgG1
Clonality	Monoclonal
MW	99998
Uniprot ID	O75473
NCBI	NP_003658.1, NP_001264155.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

