

Product Datasheet

PRMT7 antibody (orb1410050)

Description

Mouse monoclonal antibody to PRMT7

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IP, WB
Immunogen	Recombinant full-length human PRMT7 protein
Target	PRMT7
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1
Clonality	Monoclonal
MW	78.46kDa
Uniprot ID	Q9NVM4
Entrez	54496
Dilution Range	Immunoprecipitation (1-2ug per 100-500ug of total protein (1ml of cell lysate)); Western Blot (1-2ug/ml);,Flow Cytometry (1-2ug/million cells); Immunofluorescence (1-2ug/ml); ,Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.