

Product Datasheet

GTF2B/ TFIIB (Transcription Factor 2B) antibody (orb1410124)

Description

Mouse monoclonal antibody to

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IP
Immunogen	Recombinant full-length human GTF2Bprotein
Target	GTF2B
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2b
Clonality	Monoclonal
MW	34.8kDa
Uniprot ID	Q00403
Entrez	2959
Dilution Range	ELISA (For coating, order antibody without BSA); ,Flow Cytometry (1-2ug/million cells); ,Immunoprecipitation (1-2ug per 100-500ug of total protein (1ml of cell lysate)); ,Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.