

Product Datasheet

**Neuronal-Nuclei (NeuN) (Neuronal Marker)
antibody (orb1410162)**

Description

Rabbit monoclonal antibody to

Marker)

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

A synthetic peptide corresponding to residues within aa 30-60 of human NeuN protein (exact sequence is proprietary)

Target

RBFox3

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Rabbit / IgG, kappa

Clonality

Monoclonal

MW

38kDa

Uniprot ID

A6NFN3

Entrez

146713

Dilution Range

Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT), (Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95 °C followed by cooling at RT for 20 minutes), Optimal dilution for a specific application