

Product Datasheet

ACE antibody (orb1410223)

Description

Rabbit monoclonal antibody to ACE

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant fragment (around aa 652-788) of human ACE/CD143 protein
Target	ACE
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	150kDa
Uniprot ID	P12821
Entrez	1636
Dilution Range	ELISA (For coating, order antibody without BSA); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.