



## Product Datasheet

### Stathmin antibody (FITC) (orb14103)

|                            |                                                                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Stathmin antibody (FITC)                                                                                                            |
| <b>Species/Host</b>        | Rabbit                                                                                                                              |
| <b>Conjugation</b>         | FITC                                                                                                                                |
| <b>Tested Applications</b> | ICC, IF                                                                                                                             |
| <b>Immunogen</b>           | Synthetic peptide (please contact us for details)                                                                                   |
| <b>Target</b>              | Stathmin (FITC)                                                                                                                     |
| <b>Form/Appearance</b>     | Liquid: PBS, 2 mM Sodium Azide, 30% glycerol                                                                                        |
| <b>Concentration</b>       | - 100 µg (in 200 µl): 0.5 mg/ml                                                                                                     |
| <b>Storage</b>             | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.                      |
| <b>Note</b>                | For research use only                                                                                                               |
| <b>Isotype</b>             | IgG                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                          |
| <b>Purity</b>              | Polyclonal antibodies are purified by peptide affinity chromatography                                                               |
| <b>MW</b>                  | 17 kDa                                                                                                                              |
| <b>Uniprot ID</b>          | <a href="#">P54227</a> , <a href="#">P13668</a> , <a href="#">P16949</a>                                                            |
| <b>NCBI</b>                | <a href="#">NP_001138926.11</a> , <a href="#">NP_001138926.1</a> , <a href="#">NM_001145454.1</a> , <a href="#">NM_001145454.11</a> |
| <b>Entrez</b>              | <a href="#">3925</a>                                                                                                                |
| <b>Dilution Range</b>      | IF/ICC:1:50-300 (based on 0.5 mg/ml)                                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                     |



IF image of rat testis tissue using Stat...



Immunofluorescence image of rat testis t...