

Product Datasheet

TNFRSF8 (Hodgkin & Reed-Sternberg Cell Marker) antibody (orb1410338)

Description

Mouse monoclonal antibody to

Cell Marker)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

Recombinant fragment of human TNFRSF8 protein(CD30) (Epitope is cytoplasmic aa 495-595)

Target

TNFRSF8

Preservatives

200ug/ml of Ab purified by Protein A Column. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

105-120kDa

Uniprot ID

P28908

Entrez

943

Dilution Range

Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95 °C followed by cooling at RT for 20 minutes),Optimal dilution for a specific application should be determined.