

Product Datasheet

HA Tag antibody (orb1410505)

Description

Mouse monoclonal antibody to HA Tag

Species/Host

Mouse

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC-P, IP, WB

Applications
Immunogen

HA-molecule corresponding to amino acids 98-106 (YPYDVPDYA).

Target

Not Known

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8 °C. Antibody is stable for 24 months. ,Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

60kDa

Uniprot ID

Not Known

Entrez

Not Known

Dilution Range

ELISA (For coating, order antibody without BSA); Immunofluorescence (1-2ug/ml); Western Blot (1-2ug/ml); ,Immunoprecipitation (1-2ug per 100-500ug of total protein (1ml of cell lysate)); Flow Cytometry (1-2ug/million cells); ,Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95 °C followed by cooling at RT for 20 minutes),Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.