

## Product Datasheet

### AMACR antibody (orb1410512)

**Description**

Rabbit monoclonal antibody to AMACR

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

IHC-P, WB

**Applications**
**Immunogen**

Recombinant full-length human AMACR protein

**Target**

AMACR

**Preservatives**

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

**Concentration**

Purified Ab with BSA and Azide at 200ug/ml

**Storage**

Antibody with azide - store at 2 to 8 °C. Antibody is stable for 24 months. ,Non-hazardous. No MSDS required.

**Note**

For research use only

**Isotype**

Rabbit / IgG

**Clonality**

Monoclonal

**MW**

42kDa

**Uniprot ID**
**Q9UHK6**
**Entrez**
**23600**
**Dilution Range**

Western Blot (1-2ug/ml); ,Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95 °C followed by cooling at RT for 20 minutes),Optimal dilution for a specific application should be determined.

**Expiration Date**

12 months from date of receipt.