

## Product Datasheet

### BACE antibody (orb1410672)

|                            |                                        |
|----------------------------|----------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to BACE.    |
| <b>Species/Host</b>        | Rabbit                                 |
| <b>Reactivity</b>          | Human, Mouse, Rat                      |
| <b>Conjugation</b>         | Unconjugated                           |
| <b>Tested Applications</b> | WB                                     |
| <b>Note</b>                | For research use only                  |
| <b>Clonality</b>           | Polyclonal                             |
| <b>Uniprot ID</b>          | <b>P56817</b>                          |
| <b>Dilution Range</b>      | WB: WB 1:500-2000, ELISA 1:10000-20000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.        |