

Product Datasheet

MYH-pan antibody (orb1411344)

Description	Rabbit polyclonal antibody to MYH-pan.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UKX2, Q9Y623, P12882, P11055, P13533, P13535, P12883
Dilution Range	WB: Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.