

Product Datasheet

HSP70 antibody (orb1411347)

Description	Rabbit polyclonal antibody to HSP70.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P08107
Dilution Range	IF: IF: 1:50-200 Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.