

Product Datasheet

TRPV3 antibody (orb1411675)

Description	Rabbit polyclonal antibody to TRPV3.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NET8
Dilution Range	WB: Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.