

Product Datasheet

DNA pol ϵ B antibody (orb1411705)

Description	Rabbit polyclonal antibody to DNA pol ϵ B.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P56282
Dilution Range	WB: Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.