

Product Datasheet

Sox-8/9/17/18 antibody (orb1412129)

Description Rabbit polyclonal antibody to Sox-8/9/17/18.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications IHC-P, WB

Note For research use only

Clonality Polyclonal

Uniprot ID [P35713](#), [Q9H6I2](#), [P48436](#), [P57073](#)

Dilution Range WB: Western Blot: 1/500 - 1/2000.
Immunohistochemistry: 1/100 - 1/300. ELISA:
1/10000. Not yet tested in other applications.

Expiration Date 12 months from date of receipt.