

Product Datasheet

Olfactory receptor 5AP2 antibody (orb1412718)

Description	Rabbit polyclonal antibody to Olfactory receptor 5AP2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NGF4
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.