

Product Datasheet

Na⁺/K⁺-ATPase α 1 antibody (orb1412962)

Description	Rabbit polyclonal antibody to Na ⁺ /K ⁺ -ATPase α 1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P05023
Dilution Range	WB: Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.