

Product Datasheet

Trk A (Phospho-Y496) antibody (orb1414831)

Description	Rabbit polyclonal antibody to Trk A (Phospho-Y496).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04629
Dilution Range	WB: WB 1:500-2000, ELISA(peptide)1:5000-20000
Expiration Date	12 months from date of receipt.