

Product Datasheet

CaMK2 (Phospho-Thr286) antibody (orb1414833)

Description	Rabbit polyclonal antibody to CaMK2 (Phospho-Thr286).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UQM7, Q13557
Dilution Range	WB: WB 1:500-2000, ELISA(peptide)1:5000-20000
Expiration Date	12 months from date of receipt.