

Product Datasheet

Na⁺/K⁺-ATPase α 1 (Phospho-Tyr260) antibody (orb1414835)

Description	Rabbit polyclonal antibody to Na ⁺ /K ⁺ -ATPase α 1 (Phospho-Tyr260).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P05023
Dilution Range	WB: WB 1:500-2000, ELISA 1:10000-20000
Expiration Date	12 months from date of receipt.