

Product Datasheet

MDM2 (Phospho-Tyr394) antibody (orb1414851)

Description	Rabbit polyclonal antibody to MDM2 (Phospho-Tyr394).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q00987
Dilution Range	WB: WB 1:500-2000, ELISA 1:10000-20000
Expiration Date	12 months from date of receipt.