

Product Datasheet

ERK 1/2 (phospho Thr202/Y204) antibody (orb1414864)

Description	Rabbit polyclonal antibody to ERK 1/2 (phospho Thr202/Y204).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P27361, P28482
Dilution Range	IF: IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.