

Product Datasheet

NMDAε1/2 (phospho Tyr1246/1252) antibody (orb1414876)

Description	Rabbit polyclonal antibody to NMDAε1/2 (phospho Tyr1246/1252).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q13224 , Q12879
Dilution Range	IHC: Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.