

Product Datasheet

TIRAP (phospho Tyr86) antibody (orb1414884)

Description	Rabbit polyclonal antibody to TIRAP (phospho Tyr86).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P58753
Dilution Range	IHC: Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.