

Product Datasheet

Histone H2A (phospho Thr121) antibody (orb1414889)

Description	Rabbit polyclonal antibody to Histone H2A (phospho Thr121).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04908, Q7L7L0, P0C0S8, Q96KK5, P20671, Q6FI13
Dilution Range	IHC: Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.