

Product Datasheet

Na⁺/K⁺-ATPase α 1 (phospho Ser23) antibody (orb1415140)

Description	Rabbit polyclonal antibody to Na ⁺ /K ⁺ -ATPase α 1 (phospho Ser23).
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: Western Blot: 1/500 - 1/2000. IF: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.