

Product Datasheet

JIP-1 (phospho Thr103) antibody (orb1415146)

Description	Rabbit polyclonal antibody to JIP-1 (phospho Thr103).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UQF2
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.