

Product Datasheet

Smad2/3 (phospho Thr8) antibody (orb1415229)

Description	Rabbit polyclonal antibody to Smad2/3 (phospho Thr8).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P84022, Q15796
Dilution Range	WB: Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.