

## Product Datasheet

### Ephrin-B1 (phospho Tyr317) antibody (orb1415233)

|                            |                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Ephrin-B1 (phospho Tyr317).                                  |
| <b>Species/Host</b>        | Rabbit                                                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                               |
| <b>Tested Applications</b> | WB                                                                                         |
| <b>Note</b>                | For research use only                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                 |
| <b>Uniprot ID</b>          | <b>P98172</b>                                                                              |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000.<br>ELISA: 1/20000. Not yet tested in other applications. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                            |