

Product Datasheet

CaMKII α / β / δ (phospho Thr305) antibody (orb1415280)

Description	Rabbit polyclonal antibody to CaMKII α / β / δ (phospho Thr305).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UQM7 , Q13557 , Q13554
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.