

Product Datasheet

Neu (phospho Tyr1248) antibody (orb1415352)

Description	Rabbit polyclonal antibody to Neu (phospho Tyr1248).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04626
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.