

Product Datasheet

JNK1/2/3 (phospho Thr183/Y185) antibody (orb1415368)

Description	Rabbit polyclonal antibody to JNK1/2/3 (phospho Thr183/Y185).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P53779, P45983, P45984
Dilution Range	IF: IF: 1:50-200 WB 1:500-2000, IHC 1:50-300 IHC 1:50-300
Expiration Date	12 months from date of receipt.