

Product Datasheet

Histone H3 (phospho Ser28) antibody (orb1415392)

Description	Rabbit polyclonal antibody to Histone H3 (phospho Ser28).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q71DI3, P84243, P68431
Dilution Range	WB: Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.