

Product Datasheet

GSK3 β (phospho Ser9) antibody (orb1415398)

Description	Rabbit polyclonal antibody to GSK3 β (phospho Ser9).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P49841
Dilution Range	IF: IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.