

## Product Datasheet

### Phospho-Chk1 (S280) antibody (orb1415445)

|                            |                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Phospho-Chk1 (S280).                                      |
| <b>Species/Host</b>        | Rabbit                                                                                  |
| <b>Reactivity</b>          | Human                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                            |
| <b>Tested Applications</b> | WB                                                                                      |
| <b>Note</b>                | For research use only                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                              |
| <b>Uniprot ID</b>          | <b>O14757</b>                                                                           |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                         |