

Product Datasheet

α -tubulin (Acetyl Lys40) antibody (orb1415678)

Description	Mouse monoclonal antibody to α -tubulin (Acetyl Lys40).
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P68363
Dilution Range	WB: WB 1:1000-2000, IHC 1:50-100 IF 1:200
Expiration Date	12 months from date of receipt.