

Product Datasheet

P44/42 MAPK (ERK1/2) antibody (orb1415681)

Description	Mouse monoclonal antibody to P44/42 MAPK (ERK1/2).
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P27361, P28482
Dilution Range	WB: WB 1:1000-2000, IHC 1:50-100
Expiration Date	12 months from date of receipt.