

Product Datasheet

Histone H2A (Acetyl Lys15) antibody (orb1415841)

Description	Rabbit polyclonal antibody to Histone H2A (Acetyl Lys15).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	POC0S8, Q7L7L0, Q6FI13
Dilution Range	WB: WB: 1:500-1000
Expiration Date	12 months from date of receipt.